


Instruction for device driver installation

When you plug in the SensorStick for the first time, Windows will pop up a “Found New Hardware Wizard”. It will help the users easily install the device driver. We take Windows XP as an example below. Vista users will have similar procedures.


Below is a screenshot of Found New Hardware Wizard. Select “No, not this time”. Click “Next”.



Then select “Install from a list of specific location”. Click “Next”.



Unselect “search removable media” and select “include this location in the search”. Browse for the folder where you have downloaded the driver file “mchpcdc.inf”. Click “Next”.



Windows will locate the file and pop up a window showing the installation process. After it finished, the following window will show up. Just click “Finish”.



Done. Y~